

2-Way Vitrectomy Mirror Instructions



POSITION THE MIRROR

Place the 2-way mirror directly below your eyes with the slanted side that has a label pointed inwards, or towards you. Look down at the lower mirror. Your line of sight will be reflected upward into the other mirror. It will show you what is outward, and in front of you.

WHAT YOU SHOULD SEE

The mirror shows what is directly in front of you. It's not optimized to show what is higher or lower than the surface it sits on. It shows a horizontal image, so you see what is directly in front of you. Your dual sided mirror produces an image that is upright, as in normal sight.



ADJUSTING THE IMAGE

- The closer the 2-way mirror is to your face, the broader the viewing area will be.
- Tilting the 2-way mirror forward or backward will move the viewing area up or down. Although it comes pre-assembled for the best field of you possible, you can adjust the tilt of a mirror with a book, magazine, or other item. This simple, everyday angle adjustment can change the angle of the mirror and allow you to see what is higher than the surface the mirror sits on. You can also change the height or angle of the arm rest if you are using a seated support chair, if that adjusts.
- You may need to adjust the tilt of your face along with the mirror to move the image to an optimal viewing angle.
- Advanced tip: You can also get up from the seated support equipment to use the mirror to see counters and cabinets in the kitchen, but please be sure to follow your physicians instructions when spending time away from your support equipment.

